BPL_TEST2_Perfusion about:srcdoo

BPL_TEST2_Perfusion - demo

This notebook explore perfusion cultivation in comparison with ordinary continuous cultivation (chemostat) and use comparable settings to earlier notebook. Further you see here examples of interaction with the simplified commands par(), init(), simu() etc as well as direct interaction with the FMU which is called "model" here. The last simulation is always available in the workspace and called "sim_res". Note that describe() brings mainly up from descriptive information from the Modelica code from the FMU but is complemented by some information given in the Python setup file.

```
In [1]: run -i BPL_TEST2_Perfusion_explore.py
```

Linux - run FMU pre-compiled OpenModelica

Model for the process has been setup. Key commands:

- par()change of parameters and initial values
- init()change initial values only
- simu() simulate and plot
- newplot() make a new plot
- show() show plot from previous simulation
- disp()
 display parameters and initial values from the last simul
 ation
- describe() describe culture, broth, parameters, variables with value s/units


Note that both disp() and describe() takes values from the last simulation and the command process_diagram() brings up the main configuration

Brief information about a command by help(), eg help(simu) Key system information is listed with the command system_info()

```
In [2]: %matplotlib inline
   plt.rcParams['figure.figsize'] = [25/2.54, 20/2.54]
In [3]: process diagram()
```

No processDiagram.png file in the FMU, but try the file on disk.

BPL_TEST2_Perfusion about:srcdoc



```
In [4]: # Process parameters used throughout
    par(Y=0.5, qSmax=0.75, Ks=0.1)  # Cul
    par(filter_eps=0.10, filter_alpha_X=0.02, filter_alpha_S=0.10)  # Fil
    par(S_in=30.0)  # Inl
    init(V_start=1.0, VX_start=1.0)  # Pro
    eps = parDict['filter_eps']  # Pum
```

```
In [5]: # Simulation of process with flow rate close to wash-out for chemostat

init(VS_start=20)  # Process initia
    par(pump1_t1=10, pump2_t1=10)  # Pump schedule
    par(pump1_F1=2.5*0.155, pump2_F1=2.5*0.155/eps)
    par(pump1_t2=940, pump2_t2=940, pump1_t3=950, pump2_t3=950, pump1_t4=960,
    newplot(title='Perfusion cultivation - flow rate corresponding to maximal simu(60)
```

Could not find cannot import name 'dopri5' from 'assimulo.lib' (/home/janp eter/miniconda3/envs/pyfmi/lib/python3.12/site-packages/assimulo/lib/__ini t .py)

Could not find cannot import name 'rodas' from 'assimulo.lib' (/home/janpe ter/miniconda3/envs/pyfmi/lib/python3.12/site-packages/assimulo/lib/__init __.py)

Could not find cannot import name 'odassl' from 'assimulo.lib' (/home/janp eter/miniconda3/envs/pyfmi/lib/python3.12/site-packages/assimulo/lib/__ini t__.py)


Could not find ODEPACK functions.

Could not find RADAR5

Could not find GLIMDA.


BPL_TEST2_Perfusion about:srcdoc

Perfusion cultivation - flow rate corresponding to maximal rate for chemostat



In [7]: c_data=sim_res['filter.retentate.c[1]']/sim_res['filter.inlet.c[1]'] In [8]: # Simulation of process with step-wise increase of pefusion rate until wa par(pump1_F1=2.5*0.155, pump2_F1=2.5*0.155/eps) par(pump1 t2=22, pump2 t2=22)par(pump1 F2=2.5*0.35, pump2 F2=2.5*0.35/eps)par(pump1 t3=32, pump2 t3=32)par(pump1 F3=2.5*0.63, pump2 F3=2.5*0.63/eps)par(pump1_t4=42, pump2_t4=42) par(pump1 F4=2.5*0.83, pump2 F4=2.5*0.83/eps)newplot(title='Perfusion cultivation - step wise increase of perfusion ra simu(60)

BPL_TEST2_Perfusion about:srcdoc



BPL_TEST2_Perfusion

```
In [15]: model.get('filter.retentate.F')[0]
Out[15]: np.float64(-14.17499999999999)
```

Summary

- The perfusion filter had a concentration factor of cells around 1.08 and re-cycling flow was set to a factor 10 higher than the perfusion rate and changed when perfusion rate was change to keep the ratio factor 10.
- The first simulation showed that by cell retention using perfusion filter the process could be run at a perfusion flow rate at the maximal flow rate possible for corresponding chemostat culture and cell concetration increased steadily.
- The second simulation showed that with a proper startup cell concentration, the cell concentration remained constant when perfusion rate increased in a similar way as what we see in a chemostat.
- The second simulation also showed that biomass productivity in this case was increased by a factor 4.2 compared to chemostat.
- If the perfusion rate increased to higher levels washout started but the decrase of cell concentration was slow.

Some of you who read this may have your perfusion experience with CHO-cultures. For such cultures the cell concentration do increase with increase of perfusion rate and there are understood reasons for that. But for this simplified process as well as microbial processes they typically keep cell concentration constant when flow rate is chaged, and that under quite wide conditions. I will try come back to this phenomena in a later notebook.

Appendix

System information

-OS: Linux

-Python: 3.12.9

-Scipy: not installed in the notebook

-PyFMI: 2.18.0

-FMU by: OpenModelica Compiler OpenModelica 1.26.0~dev-200-gcb3254b

-FMI: 2.0

-Type: FMUModelME2

-Name: BPL.Examples_TEST2.Perfusion -Generated: 2025-07-28T07:59:36Z

-MSL: 4.1.0

-Description: Bioprocess Library version 2.3.1

-Interaction: FMU-explore version 1.0.0

In [21]: !lsb_release -a

No LSB modules are available.

Distributor ID: Ubuntu

Description: Ubuntu 24.04.3 LTS

Release: 24.04 Codename: noble