BPL_TEST2_Perfusion script with FMPy

The key library FMPy is installed.

After the installation a small application BPL_TEST2_Perfusion is loaded and run. You can continue with this example if you like.

```
!lsb_release -a # Actual VM Ubuntu version used by Google
No LSB modules are available.
    Distributor ID: Ubuntu
                    Ubuntu 22.04.3 LTS
    Description:
                    22.04
    Release:
    Codename:
                    jammy
%env PYTHONPATH=
→ env: PYTH0NPATH=
!wget https://repo.anaconda.com/miniconda/Miniconda3-py312_24.3.0-0-Linux-x86_64.sh
!chmod +x Miniconda3-py312_24.3.0-0-Linux-x86_64.sh
!bash ./Miniconda3-py312_24.3.0-0-Linux-x86_64.sh -b -f -p /usr/local
sys.path.append('/usr/local/lib/python3.12/site-packages/')
   --2024-05-15 08:14:37-- https://repo.anaconda.com/miniconda/Miniconda3-py312 24.3.0-0-Linux-
    Resolving repo.anaconda.com (repo.anaconda.com)... 104.16.191.158, 104.16.32.241, 2606:4700::
    Connecting to repo.anaconda.com (repo.anaconda.com)|104.16.191.158|:443... connected.
    HTTP request sent, awaiting response... 200 OK
    Length: 143351488 (137M) [application/octet-stream]
    Saving to: 'Miniconda3-py312_24.3.0-0-Linux-x86_64.sh'
    Miniconda3-py312_24 100%[===========] 136.71M 89.5MB/s
                                                                         in 1.5s
    2024-05-15 08:14:39 (89.5 MB/s) - 'Miniconda3-py312_24.3.0-0-Linux-x86_64.sh' saved [14335148]
    PREFIX=/usr/local
    Unpacking payload ...
    Installing base environment...
    Preparing transaction: ...working... done
    Executing transaction: ...working... done
    installation finished.
!conda update -n base -c defaults conda --yes
→ Channels:
     defaults
    Platform: linux-64
    Collecting package metadata (repodata.json): done
    Solving environment: done
    ## Package Plan ##
      environment location: /usr/local
      added / updated specs:

    conda
```

The following packages will be downloaded:

package	build	
conda-24.5.0 frozendict-2.4.2 openssl-3.0.13	py312h06a4308_0 py312h06a4308_0 h7f8727e_1	1.2 MB 36 KB 5.2 MB
	Total:	6.5 MB

The following NEW packages will be INSTALLED:

frozendict pkgs/main/linux-64::frozendict-2.4.2-py312h06a4308_0

The following packages will be UPDATED:

conda 24.3.0-py312h06a4308_0 --> 24.5.0-py312h06a4308_0 openssl 3.0.13-h7f8727e_0 --> 3.0.13-h7f8727e_1

Downloading and Extracting Packages: openssl-3.0.13 | 5.2 MB | : 0% 0/1 [00:00<?, ?it/s] conda-24.5.0 | 1.2 MB 0% 0/1 [00:00<?, ?it/s] | : frozendict-2.4.2 | 36 KB 0% 0/1 [00:00<?, ?it/s] openssl-3.0.13 | 5.2 MB 0% 0.002997347135570501/1 [00:00<00:50, 50.71s/it] | : frozendict-2.4.2 | 36 KB | : 44% 0.43853215920344746/1 [00:00<00:00, 2.77it/s] frozendict-2.4.2 | : 100% 1.0/1 [00:00<00:00, 2.77it/s] | 36 KB conda-24.5.0 | 1.2 MB | : 100% 1.0/1 [00:00<00:00, 1.84it/s]

Preparing transaction: done Verifying transaction: done Executing transaction: done

!conda --version
!python --version

conda 24.5.0 Python 3.12.2

!conda install -c conda-forge fmpy --yes # Install the key package

 $\overline{\mathbf{T}}$

```
Preparing transaction: done
Verifying transaction: done
Executing transaction: done

#!conda install -c conda-forge matplotlib --yes

#!conda install -c conda-forge scipy --yes

#!conda install -c conda-forge openpyxl --yes

#!conda install -c conda-forge xlrd --yes
```

Notes of BPL_TEST2_Perfusion

This notebook explore perfusion cultivation in comparison with ordinary continuous cultivation (chemostat) and use comparable settings to earlier notebook. Further you see here examples of interaction with the simplified commands par(), init(), simu() etc as well as direct interaction with the FMU which is called "model" here. The last simulation is always available in the workspace and called "sim_res". Note that describe() brings mainly up from descriptive information from the Modelica code from the FMU but is complemented by some information given in the Python setup file.

Now specific installation run a simulation and notebook for that Start with connecting to Github. Then upload the two files:

- FMU BPL_TEST2_Perfusion_linux_om_me.fmu
- Setup-file BPL_TEST2_Perfusion_fmpy_explore.py

%%bash

git clone https://github.com/janpeter19/BPL_TEST2_Perfusion

Cloning into 'BPL_TEST2_Perfusion'...

%cd BPL_TEST2_Perfusion

/content/BPL_TEST2_Perfusion

run -i BPL_TEST2_Perfusion_fmpy_explore.py

→ Linux - run FMU pre-comiled OpenModelica 1.21.0

Model for bioreactor has been setup. Key commands:

- par() change of parameters and initial values
- init() change initial values only
- simu() simulate and plot
- newplot() make a new plot
 show() show plot from previous
- show()- show plot from previous simulation- disp()- display parameters and initial values from the last simulation
- describe() describe culture, broth, parameters, variables with values/units


Note that both disp() and describe() takes values from the last simulation and the command process_diagram() brings up the main configuration

Brief information about a command by help(), eg help(simu)
Key system information is listed with the command system_info()

%matplotlib inline
plt.rcParams['figure.figsize'] = [25/2.54, 20/2.54]

process_diagram()


No processDiagram.png file in the FMU, but try the file on disk.




```
# Process parameters used throughout
par(Y=0.5, qSmax=0.75, Ks=0.1)  # Culture
par(filter_eps=0.10, filter_alpha_X=0.02, filter_alpha_S=0.10)  # Filter
par(S_in=30.0)  # Inlet substrate conc
init(V_start=1.0, VX_start=1.0)  # Process initial conditions
eps = parDict['filter_eps']  # Pump schedule parameter
```

Simulation of process with flow rate clot to wash-out for chemostat

newplot(title='Perfusion cultivation - flow rate corresponding to maximal rate for chemostat')
simu(60)




Simulation of process with flow rate close to wash-out for chemostat



```
# Concentration factor of the filter
c=model_get('filter.retentate.c[1]')/model_get('filter.inlet.c[1]')
print('Conc factor of perfusion filter =', np.round(c,3))
Example 2.186 Conc factor of perfusion filter = 1.186
c_data=sim_res['filter.retentate.c[1]']/sim_res['filter.inlet.c[1]']
print('Conc factor variation', np.round(min(c_data[151:]), 3), np.round(max(c_data[151:]),3))

→ Conc factor variation 1.186 1.186
# Simulation of process with step-wise increase of pefusion rate until wash-out.
# This means that re-circulation rate change at the same time as the perfusion rate.
init(VS_start=150)
                                                          # Process initial varied
par(pump1_t1=12, pump2_t1=12)
                                                          # Pump schedule - recycle flow 10 times pe
par(pump1_F1=2.5*0.155, pump2_F1=2.5*0.155/eps)
par(pump1_t2=22, pump2_t2=22)
par(pump1_F2=2.5*0.35, pump2_F2=2.5*0.35/eps)
par(pump1_t3=32, pump2_t3=32)
par(pump1_F3=2.5*0.63, pump2_F3=2.5*0.63/eps)
par(pump1_t4=42, pump2_t4=42)
par(pump1 F4=2.5*0.83, pump2 F4=2.5*0.83/eps)
```

newplot(title='Perfusion cultivation - step wise increase of perfusion rate and max productivity 4>
simu(60)



```
# Simulation without a plot and just to check typical values at high production rate #simu(40)

#c_data=sim_res['filter.retentate.c[1]']/sim_res['filter.inlet.c[1]']

#print('Conc factor variation', np.round(min(c_data[190:]), 3), 'to', np.round(max(c_data[190:]),

#describe('cstrProdMax')

# The maximal biomass productivity before washout is obtained aroudn 40 hours np.round(model_get('harvesttank.inlet.F')*model_get('harvesttank.inlet.c[1]'),1)

$\sum_{\text{2}} 25.2

# Thus perfusion (with this filter) brings a productivity improvement of about np.round(23.5/5.6,1)

$\sum_{\text{2}} 4.2

# Finally we check the filter flow rates at time 40 hour - note the negative sign for outlfow model_get('filter.inlet.F')

$\sum_{\text{2}} 20.749999999999996
```

Summary

- The perfusion filter had a concentration factor of cells around 1.08 and re-cycling flow was set to a factor 10 higher than the perfusion rate and changed when perfusion rate was change to keep the ratio factor 10.
- The first simulation showed that by cell retention using perfusion filter the process could be run at a perfusion flow rate at the maximal flow rate possible for corresponding chemostat culture and cell concetration increased steadily.
- The second simulation showed that with a proper startup cell concentration, the cell concentration remained constant when perfusion rate increased in a similar way as what we see in a chemostat.
- The second simulation also showed that biomass productivity in this case was increased by a factor 4.2 compared to chemostat.
- If the perfusion rate increased to higher levels washout started but the decrase of cell concentration was slow.

Some of you who read this may have your perfusion experience with CHO-cultures. For such cultures the cell concentration do increase with increase of perfusion rate and there are understood reasons for that. But for this simplified process as well as microbial processes they typically keep cell concentration constant when flow rate is chaged, and that under quite wide conditions. I will try come back to this phenomena in a later notebook.

```
# List of components in the process setup and also a couple of other things like liquidphase and D
describe('parts')
['bioreactor', 'bioreactor.culture', 'D', 'feedtank', 'filter', 'harvesttank', 'schemePump1',
describe('MSL')
🚁 MSL: 3.2.3 - used components: RealInput, RealOutput, CombiTimeTable, Types
system_info()
    System information
     -OS: Linux
     -Python: 3.10.12
     -Scipy: not installed in the notebook
     -FMPy: 0.3.20
     -FMU by: OpenModelica Compiler OpenModelica 1.21.0
     -FMI: 2.0
     -Type: ME
     -Name: BPL_TEST2.Perfusion
     -Generated: 2024-03-05T09:10:09Z
     -MSL: 3.2.3
     -Description: Bioprocess Library version 2.2.0
     -Interaction: FMU-explore for FMPy version 1.0.0
```